


In-Ovation L Patient Benefits


- Absolutely undetectable
- Small, comfortable to wear
- Speak clearly and confidently
- Extended intervals between appointments
- Short treatment time
- Minimal chair time



- 1 True twin design
- 2 Triple Chamfered slot walls
- 3 Smooth swept tie-wings reduce occlusal interference
- 4 Interactive spring clip provides full slot coverage with two point contact for true twin rotational control and is gold for enhanced visibility
- 5 Anatomically correct base



In-Ovation's Interactive Control




Experience the difference of with In-Ovation L

Dentsply GAC has translated the proven clinical benefits of In-Ovation into a lingual design. Only an interactive bracket provides full functionality throughout the course of treatment. This allows you to enjoy the early benefits of a free-sliding, passive design and the precise control of an active one in the later stages of treatment. The In-Ovation L self-ligating bracket gives you a complete range of control options... simply by changing the archwire.


Lingual Mushroom & Universal Arch Form Matrix

Sentalloy/Neo-Sentalloy Maxillary



Item#	Description	Force	Arch Form	Item#	Description	Force	Arch Form	Item#	Description	Force	Arch Form	Item#	Description	Force	Arch Form
02-911-101	SENTALLOY .012	MED	SM	02-911-102	SENTALLOY .012	MED	MED	02-911-103	SENTALLOY .012	MED	LG	02-911-104	SENTALLOY .012	MED	XLG
02-911-111	SENTALLOY .014	MED	SM	02-911-112	SENTALLOY .014	MED	MED	02-911-113	SENTALLOY .014	MED	LG	02-911-114	SENTALLOY .014	MED	XLG
02-923-131	NEO SENT .017X.025	MED	SM	02-923-132	NEO SENT .017X.025	MED	MED	02-923-133	NEO SENT .017X.025	MED	LG	02-923-134	NEO SENT .017X.025	MED	XLG


Mandibular



Item#	Description	Force	Arch Form	Item#	Description	Force	Arch Form	Item#	Description	Force	Arch Form	Item#	Description	Force	Arch Form
02-911-601	SENTALLOY .012	MED	SM	02-911-602	SENTALLOY .012	MED	MED	02-911-603	SENTALLOY .012	MED	LG	02-911-604	SENTALLOY .012	MED	XLG
02-911-611	SENTALLOY .014	MED	SM	02-911-612	SENTALLOY .014	MED	MED	02-911-613	SENTALLOY .014	MED	LG	02-911-614	SENTALLOY .014	MED	XLG
02-923-631	NEO SENT .017X.025	MED	SM	02-923-632	NEO SENT .017X.025	MED	MED	02-923-633	NEO SENT .017X.025	MED	LG	02-923-634	NEO SENT .017X.025	MED	XLG

RESOLVE[®] Beta-Titanium Universal Arch Wires

03-016-000 RESOLVE .016 LING UNIV ARCH
03-622-000 RESOLVE .016X022 LING UNIV ARCH



The Perfect Aesthetic Solution

More and more, adults want to experience the benefits of a perfect smile, without the obvious appearance of metal or ceramic brackets. Now you have a completely undetectable solution to offer your patients with In-Ovation[®] L.



Unprecedented Ease of Use

Lingual brackets have been around for 30 years, dramatically changing the concept of aesthetic brackets. In-Ovation L leverages all the advantages of the In-Ovation self-ligating line of brackets to solve one of the biggest problems doctors have traditionally encountered with lingual systems — archwire ligation. The In-Ovation L clip opens and closes easily and quickly to make archwire ligation literally a “snap.”

Doctor and Staff Benefits:

Feature

- Low profile
- Beveled slot walls
- Interactive Clip
- True Twin Design
- Anatomically correct Base
- Metal Injection Molded (MIM)
- CNC milled

Benefits


- Patient comfort and avoiding gingival interference for better patient experience
- Facilitates wire engagement and superior sliding mechanics
- Provides full functionality throughout the course of treatment for greater control
- Offers optimal mesial, distal span for achieving superior rotation control
- Design Adapts smoothly to the lingual anatomy
- Small size and greater strength
- True slot dimensions and smooth finish



Part Art. Part Science. All Orthodontics.

Enhancing the Art and Science of Orthodontics

From esthetics perfection to essential care, Orthodontics is about moving the individual parts to perfect the collective whole. It's why DENTSPLY GAC is moving the field of orthodontics forward to deliver the highest in patient care. From brackets, force-intuitive wires and highly efficient buccal tubes to bonding, auxiliary appliances and innovative software programs, our mission is to help you achieve rapid, optimal treatment results while streamlining your orthodontic practice. To ask your GAC representative for more information, please call 800.645.5530, or visit www.gacintl.com.



Part Art. Part Science. All Orthodontics.

800.645.5530 | www.gacintl.com



Leading the Way in Self-Ligating Lingual Brackets

800.645.5530 | www.gacintl.com

© 2011 Dentsply GAC International. All Rights Reserved. 2.2011 • U.S. Patent: 6,776,613. Other U.S. and foreign patents pending and/or approved. 120-090-01

For better dentistry
DENTSPLY
GAC

For better dentistry
DENTSPLY
GAC