

From Great Resolve Comes Great Results

The name is no accident. Resolve is more than a product line. It's a mindset that permeates everything we do—from the frontline of product innovation to the bottom line of every orthodontist we serve. Let's take a moment to consider the origins.

Resolve (noun)

1. Firmness of purpose; resolution.
2. A determination or decision; a fixed purpose.
3. A formal resolution made by a deliberative body.
4. Boldness, courage, fixed purpose, intention, objective, project, purpose, purposefulness, resolution, steadfastness, will, willpower.

Resolve

Resolve to SUCCEED. Resolve to SERVE our customers. Resolve to ADVANCE orthodontics.

From the newly introduced Accu-Force™ and Shape-Fit™ products to the new Resolve® line of brackets, these products are the result of considerable dedication, determination and the resolve to overcome a challenge like no other we've faced. The Resolve name represents both the culmination of a goal to prevail and a starting point for a continual journey. And while the Resolve line is but one of the many products that make up our complete offering, we feel it is an accurate representation of who we are and what we stand for.

Part Art. Part Science. All Resolve.



 **resolve**[®]
VS Bracket

- Low profile, mini bracket designed with patient comfort in mind
- Vertical slot feature for use of auxiliaries
- Compound contoured base provides an anatomical fit
- Metal Injection Molding (MIM) design

2011 ©DENTSPLY GAC Intl. All Rights Reserved. 120-172-02 2011.6


Rx Only

DENTSPLY INTERNATIONAL
DENTSPLY GAC INTERNATIONAL
355 Knickerbocker Avenue
Bohemia, NY 11716
1-888-422-4685
Outside US: 1-631-419-1700
www.gacintl.com

Part Art. Part Science. All Orthodontics.

For better dentistry
DENTSPLY
GAC

Strength and Control with a Mini-Twin Bracket System



System Benefits

- Low profile bracket with rounded edges and a polished finish are designed with patient comfort and esthetics in mind
- Ample tie-wings allow for easy ligation
- Compound contour base provides an anatomical fit
- Available in Roth, Cetlin, and Bi-Dimensional prescriptions

Vertical Slot

The vertical slot design accommodates the use of auxiliaries during treatment.

Micro-Etched Base

The Resolve VS bracket uses an 80 gauge, single-mesh micro-etched base, designed for predictable bonding.



AVAILABLE PRESCRIPTIONS

Roth Rx

Maxillary Arch				.018				.022			
Tooth	T	A	OF	Right	ID	Left	ID	Right	ID	Left	ID
U1	12T	5A	00	172-111-03	●	172-211-03	●	172-112-03	●	172-212-03	●
U2	8T	9A	00	172-121-03	●	172-221-03	●	172-122-03	●	172-222-03	●
U3/distal hook	-2T	7A	00	172-131-13	●	172-231-13	●	172-132-13	●	172-232-13	●
U4-5	-7T	0A	00	172-141-03	●	172-241-03	●	172-142-03	●	172-242-03	●
Mandibular Arch				.018				.022			
L1-2	0T	0A	00	172-311-03	●	172-411-03	●	172-312-03	●	172-412-03	●
L3/distal hook	-11T	5A	00	172-331-10	●	172-431-10	●	172-332-10	●	172-432-10	●
L4-5	-17T	0A	00	172-341-00	○	172-441-00	○	172-342-00	○	172-442-00	○
L5	-22T	0A	00	172-351-00	●	172-451-00	●	172-352-00	●	172-452-00	●

Cetlin Rx

Maxillary Arch				.018				.022			
Tooth	T	A	OF	Right	ID	Left	ID	Right	ID	Left	ID
U1	12T	5A	00	172-111-03	●	172-211-03	●	172-112-03	●	172-212-03	●
U2	8T	9A	00	172-121-03	●	172-221-03	●	172-122-03	●	172-222-03	●
U3/distal hook	-2T	7A	00	172-131-13	●	172-231-13	●	172-132-13	●	172-232-13	●
U4-5	-7T	0A	00	172-141-03	●	172-241-03	●	172-142-03	●	172-242-03	●
Mandibular Arch				.018				.022			
L1-2	0T	0A	00	172-311-03	●	172-411-03	●	172-312-03	●	172-412-03	●
L3/distal hook	-5T	7A	00	172-331-13	●	172-431-13	●	172-332-13	●	172-432-13	●
L4-5	-11T	0A	00	172-341-03	○	172-441-03	○	172-342-03	○	172-442-03	○
L4-5	-17T	0A	00	172-341-00	○	172-341-00	○	172-342-00	○	172-342-00	○

Bi-Dimensional Rx

Maxillary Arch				.018				.022			
Tooth	T	A	OF	Right	ID	Left	ID	Right	ID	Left	ID
U1	12T	5A	00	172-111-03	●	172-211-03	●				
U2	8T	9A	00	172-121-03	●	172-221-03	●				
U3/distal hook	0T	7A	00					172-132-14	●	172-232-14	●
U4-5	0T	0A	00					172-142-04	●	172-242-04	●
Mandibular Arch				.018				.022			
L1-2	0T	0A	00	172-311-03	●	172-411-03	●				
L3/distal hook	0T	7A	00					172-332-14	●	172-332-14	●
L4-5	0T	0A	00					172-342-04	●	172-342-04	●