

Bands

Selecting the proper orthodontic band is crucial in the corrective process. The criteria for selecting the proper band can include ease and accuracy of placement, elasticity, and how the band contours to the anatomy of the tooth. In addition to a brilliant, hand-polished finish, our bands have either acid-etched or micro-etched

properties that enhance retention to the tooth. Our bands also have a convenient labeling system that is laser-etched right onto the surface. At DENTSPLY GAC we understand how crucial the bands are in the anchorage function, and we offer a comprehensive line of orthodontic brackets that are guided by these criteria.

For better dentistry

DENTSPLY
GAC

800.645.5530 | www.gacintl.com

Bands Section Content

Band Identification	77
Snap-Fit Molar Bands	78
Ideal Molar Bands	79
Ovation Molar Bands	80
Snap-Fit Bicuspid Bands	81
Snap-Fit Cuspid Bands	82
Snap-Fit Anterior Bands	82
Band Material & Instruments	83

Band Identification

DENTSPLY GAC offers three categories of Bands – Snap-Fit, Ideal Molar, and Ovation. Please refer the guide below for band identification:

Snap-Fit Bands

First Molar Example

U L 1

Coding Information

“U” indicates upper, “L” indicates lower
“L” indicates left, “R” indicates right
“1” indicates size, sizes are available from 1 through 32



Second Molar Example

S U L 1

Coding Information

“S” indicates second molar
“U” indicates upper, “L” indicates lower
“L” indicates left, “R” indicates right
“1” indicates size, sizes are available from 1 through 30



Ideal Molar Bands

First Molar Example

5 U L

Coding Information

“5” indicates size, sizes are available from 5 through 38
“U” indicates upper, “L” indicates lower
“L” indicates left, “R” indicates right



Second Molar Example

1 S U L

Coding Information

“1” indicates size, sizes are available from 1 through 32
“S” indicates second molar
“U” indicates upper, “L” indicates lower
“L” indicates left, “R” indicates right



Ovation Bands

First Molar Example

U 6 L 1

Coding Information

“U” indicates upper, “L” indicates lower
“6” indicates 1st molar
“L” indicates left, “R” indicates right
“1” indicates size, sizes are available from 1 through 32



Snap-Fit Bicuspid Bands

Upper Example

U 1

Coding Information

“U” indicates upper
“1” indicates size, sizes are available from 1 through 32



Lower First Example

L 1

Coding Information

“L” indicates lower
“1” indicates size, sizes are available from 1 through 25



Lower Second Example

32 L

Coding Information

“32” indicates size, sizes are available from 6 through 32
“L” indicates lower



Snap-Fit Cuspid Bands

Cuspid Example

L 1

Coding Information

“U” indicates upper, “L” indicates lower
“1” indicates size, sizes are available from 000 through 16



Snap-Fit Anterior Bands

Cuspid Example

U 9

Coding Information

“U” indicates universal
“9” indicates size, sizes are available from 0 through 24



Bands

Snap-Fit Molar Bands

Our Snap-Fit Molar Bands offer greater elasticity and superior fit. GAC's exclusive process produces bands with improved morphology, elasticity, and memory. After fitting and removal, the band returns to its original size, assuring that the cemented band is as accurate as the initial fit. Snap-Fit Molar Bands are offered in graduated, proportional sizes, maxillary and mandibular, left and right (32 sizes for first molars, 30 sizes for second molars). Plain or pre-welded assemblies, including our full line of Micro MIM buccal tubes, are offered individually and in convenient kits. All Snap-Fit Bands are marked using our patented Electro Chemical Etching process for easy-to-read permanent identification, and interiors are acid-etched for roughness to increase bond strength.



Snap-Fit First Molar Distribution Kit - Distribution by Quadrant: (UR,UL,LR,LL)

Band Size	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18	19	20	21	22	23	24	25	26	27	28	29	30	31	32
Starter											1	2	2	2	2	2	2	2	2	2	2	2	2									
Intro	1	1	1	1	1	1	1	1	1	1	2	2	2	2	2	3	3	3	3	2	2	2	2	2	1	1	1	1	1	1	1	1
Standard	1	1	1	1	2	2	2	3	4	4	4	5	5	6	6	6	6	6	6	5	5	4	3	2	2	2	1	1	1	1	1	
Master	1	1	1	2	2	2	4	4	6	6	8	10	10	12	14	16	16	16	14	12	10	8	6	4	4	2	2	2	2	1	1	1

Snap-Fit Second Molar Distribution Kit - Distribution by Quadrant: (UR,UL,LR,LL)

Band Size	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18	19	20	21	22	23	24	25	26	27	28	29	30	31	32
Starter											1	2	2	2	2	2	2	2	2	2	2	2	2									
Intro	1	1	1	1	1	1	1	1	1	2	2	2	2	2	2	3	3	3	3	2	2	2	2	2	2	1	1	1	1	1		
Standard	1	1	1	2	2	2	2	3	4	4	4	5	5	6	6	6	6	6	6	5	5	4	3	2	2	2	2	1	1	1		
Master	1	1	1	2	2	2	4	4	6	6	8	10	10	12	14	16	16	16	16	12	10	8	6	4	4	2	2	2	2	1		



Upper First Molar, Plain

Reference Number

Starter Kit: 50 Assorted Bands, 13 Sizes	KIT26-050-16
Intro Kit: 100 Assorted Bands, 32 Sizes	KIT26-100-16
Standard Kit: 200 Assorted Bands, 32 Sizes	KIT26-200-16
Master Kit: 400 Assorted Bands, 32 Sizes	KIT26-400-16



Lower First Molar, Plain

Reference Number

Starter Kit: 50 Assorted Bands, 13 Sizes	KIT26-050-26
Intro Kit: 100 Assorted Bands, 32 Sizes	KIT26-100-26
Standard Kit: 200 Assorted Bands, 32 Sizes	KIT26-200-26
Master Kit: 400 Assorted Bands, 32 Sizes	KIT26-400-26



Upper Second Molar, Plain

Reference Number

Starter Kit: 50 Assorted Bands, 13 Sizes	KIT26-050-17
Intro Kit: 100 Assorted Bands, 30 Sizes	KIT26-100-17
Standard Kit: 200 Assorted Bands, 30 Sizes	KIT26-200-17
Master Kit: 400 Assorted Bands, 30 Sizes	KIT26-400-17



Lower Second Molar, Plain

Reference Number

Starter Kit: 50 Assorted Bands, 13 Sizes	KIT26-050-27
Intro Kit: 100 Assorted Bands, 30 Sizes	KIT26-100-27
Standard Kit: 200 Assorted Bands, 30 Sizes	KIT26-200-27
Master Kit: 400 Assorted Bands, 30 Sizes	KIT26-400-27

See the buccal tube section of this catalog for Micro MIM attachments for your technique.

Ideal Molar Bands

Ideal Molar Bands adapt to the contour of the tooth as they are positioned. These bands are easier to re-contour than our Snap-Fit Bands. In addition, users of our Ideal Bands will enjoy an edge hardness that prevents edge crushing, as well as a crown height. They are made from seamless stainless steel and are stretched minimally during contouring to assure greater strength. The interiors of all bands are etched for roughness to increase bond strength. Sized 5-38 for first molars and 1-32 for second molars, these “softer” bands provide excellent fit and stability. They are anatomically contoured for left and right, and they are marked using our patented Electro Chemical Etch process for easy-to-read permanent identification, and the interiors are acid-etched for roughness to increase bond strength. This process provides permanent band identification that withstands sterilization and will not weaken the integrity of the band nor stain nor anneal, as laser identification can.



Ideal First Molar Distribution Kit - Distribution by Quadrant: (UR,UL,LR,LL)

Band Size	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18	19	20	21	22	23	24	25	26	27	28	29	30	31	32	33	34	35	36	37	38	
Starter												1	2	2	2	2	2	2	2	2	2	2	2	2	2	2													
Intro				1	1	1	1	1	1	1	1	1	1	2	2	2	2	2	3	3	3	2	2	2	2	2	1	1	1	1	1	1	1	1	1	1	1	1	1
Standard				1	1	1	1	2	2	3	4	4	4	5	6	6	6	6	6	6	5	5	4	4	3	3	2	1	1	1	1	1	1	1	1	1	1	1	1
Master				2	2	2	2	4	4	6	8	8	8	10	12	12	12	12	12	12	10	10	8	8	6	6	4	2	2	2	2	2	2	2	2	2	2	2	2

Ideal Second Molar Distribution Kit - Distribution by Quadrant: (UR,UL,LR,LL)

Band Size	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18	19	20	21	22	23	24	25	26	27	28	29	30	31	32						
Starter											1	2	2	2	2	2	2	2	2	2	2	2	2															
Intro	1	1	1	1	1	1	1	1	1	1	1	1	2	2	2	2	2	2	3	3	3	3	2	2	2	2	1	1	1	1	1	1	1	1	1	1	1	1
Standard	1	1	1	1	2	2	2	3	4	4	4	5	5	6	6	6	6	6	6	5	5	4	3	2	2	2	1	1	1	1	1	1	1	1	1	1	1	1
Master	1	1	1	2	2	2	4	4	6	6	8	10	10	12	14	16	16	16	14	12	10	8	6	4	4	2	2	2	2	2	1	1	1	1	1	1	1	1

Upper First Molar, Plain

Reference Number

Starter Kit: 50 Assorted Bands, 13 Sizes	KIT16-050-11
Intro Kit: 100 Assorted Bands, 34 Sizes	KIT16-100-11
Standard Kit: 200 Assorted Bands, 34 Sizes	KIT16-200-11
Master Kit: 400 Assorted Bands, 34 Sizes	KIT16-400-11



Lower First Molar, Plain

Reference Number

Starter Kit: 50 Assorted Bands, 13 Sizes	KIT16-050-21
Intro Kit: 100 Assorted Bands, 34 Sizes	KIT16-100-21
Standard Kit: 200 Assorted Bands, 34 Sizes	KIT16-200-21
Master Kit: 400 Assorted Bands, 34 Sizes	KIT16-400-21



Upper Second Molar, Plain

Reference Number

Starter Kit: 50 Assorted Bands, 13 Sizes	KIT16-050-17
Intro Kit: 100 Assorted Bands, 32 Sizes	KIT16-100-17
Standard Kit: 200 Assorted Bands, 32 Sizes	KIT16-200-17
Master Kit: 400 Assorted Bands, 32 Sizes	KIT16-400-17



Lower Second Molar, Plain

Reference Number

Starter Kit: 50 Assorted Bands, 13 Sizes	KIT16-050-27
Intro Kit: 100 Assorted Bands, 32 Sizes	KIT16-100-27
Standard Kit: 200 Assorted Bands, 32 Sizes	KIT16-200-27
Master Kit: 400 Assorted Bands, 32 Sizes	KIT16-400-27



Bands

Ovation Molar Bands

Ovation bands were designed with a shorter crown height and evaluated to satisfy the needs of the Roth Williams International Society of Orthodontists. The soft contour and easier sizing make placement and adaptation simple. The strength is similar to previous bands designed for the Roth Williams International Society, and they hold up well throughout treatment. Ovation bands are available in first molar sizes 1 through 32. Also GAC's triple-pulse welding process fuses tubes to the bands, and the interiors are micro-etched for enhanced retention.



Ovation First Molar Distribution Kit - Distribution by Quadrant: (UR,UL,LR,LL)

Band Size	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18	19	20	21	22	23	24	25	26	27	28	29	30	31	32
Starter											1	2	2	2	2	2	2	2	2	2	2	2	2									
Intro	1	1	1	1	1	1	1	1	1	1	2	2	2	2	2	3	3	3	3	2	2	2	2	2	1	1	1	1	1	1	1	1
Standard	1	1	1	1	2	2	2	3	4	4	4	5	5	6	6	6	6	6	6	5	5	4	3	2	2	2	1	1	1	1	1	
Master	1	1	1	2	2	2	4	4	6	6	8	10	10	12	14	16	16	16	14	12	10	8	6	4	4	2	2	2	2	1	1	1



Upper First Molar, Plain

Reference Number

Starter Kit: 50 Assorted Bands, 13 Sizes

KIT36-050-16

Intro Kit: 100 Assorted Bands, 32 Sizes

KIT36-100-16

Standard Kit: 200 Assorted Bands, 32 Sizes

KIT36-200-16

Master Kit: 400 Assorted Bands, 32 Sizes

KIT36-400-16



Lower First Molar, Plain

Reference Number

Starter Kit: 50 Assorted Bands, 13 Sizes

KIT36-050-26

Intro Kit: 100 Assorted Bands, 32 Sizes

KIT36-100-26

Standard Kit: 200 Assorted Bands, 32 Sizes

KIT36-200-26

Master Kit: 400 Assorted Bands, 32 Sizes

KIT36-400-26

Snap-Fit Bicuspid Bands

Snap-Fit Bicuspid Bands complement our molar bands with superior fit and morphology. In upper and lower, our bicuspid bands are carefully proportioned in height and circumference. Their improved morphology and snap-fit design ensure easier seating. Bands are seated into position with a rocking motion using alternating mesial/distal pressure. Our easy-to-fit, seamless steel bands are offered in 32 graduated sizes of upper bicuspid. There are 25 sizes in lower first and 27 sizes in lower second bicuspid bands. Plain or with pre-welded assemblies, bands are available individually or in economical professional kits. All bands are marked using our patented Electro Chemical Etching process for easy-to-read permanent identification, and the interiors are acid-etched for roughness to increase bond strength.



Snap-Fit Upper Bicuspid Band Kit Distribution

Band Size	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18	19	20	21	22	23	24	25	26	27	28	29	30	31	32
Starter						2	2	2	2	2	2	2	2	2	2	4	4	4	2	2	2	2	2	2	2	2						
Intro	2	2	2	2	2	2	2	2	2	2	4	4	4	4	6	6	6	6	6	4	4	4	4	2	2	2	2	2	2	2	2	2
Standard	2	2	2	2	2	4	6	6	6	6	6	8	8	10	12	14	14	14	12	10	8	8	6	6	6	4	4	4	2	2	2	2
Master	4	4	4	6	6	10	10	14	14	16	16	16	18	18	20	22	22	22	20	18	16	16	14	14	12	12	8	8	8	4	4	4

Snap-Fit Lower First Bicuspid Band Kit Distribution

Band Size	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18	19	20	21	22	23	24	25
Starter			2	4	4	4	4	4	4	4	4	4	4	4	4										
Intro	2	2	2	2	2	2	4	4	4	4	4	6	6	8	8	8	8	4	4	4	4	2	2	2	2
Standard	2	2	2	2	4	4	8	8	8	8	10	12	12	18	18	18	18	10	8	8	8	4	4	2	2
Master	4	4	4	6	6	8	8	14	14	18	18	24	24	36	36	36	36	24	24	16	16	8	8	4	4

Snap-Fit Lower Second Bicuspid Band Kit Distribution

Band Size	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18	19	20	21	22	23	24	25	26	27	28	29	30	31	32
Starter						1	1	2	2	2	2	2	2	2	2	2	2	2	2	2	2	2	2	2	2	2	2	2	2	2	1	1
Intro						2	2	2	2	2	2	2	4	4	4	6	6	8	8	8	6	6	4	4	4	2	2	2	2	2	2	2
Standard						2	2	2	4	4	6	6	6	8	10	12	14	16	16	16	14	12	10	8	6	6	6	4	4	2	2	2
Master						2	4	4	6	8	10	12	12	16	20	26	28	34	34	34	28	26	20	16	14	12	10	8	6	4	4	2

Upper Bicuspid, Plain

Starter Kit: 50 Assorted Bands, 22 Sizes

Intro Kit: 100 Assorted Bands, 32 Sizes

Standard Kit: 200 Assorted Bands, 32 Sizes

Master Kit: 400 Assorted Bands, 32 Sizes

Reference Number

KIT15-050-11

KIT15-100-11

KIT15-200-11

KIT15-400-11



Lower First Bicuspid, Plain

Starter Kit: 50 Assorted Bands, 13 Sizes

Intro Kit: 100 Assorted Bands, 25 Sizes

Standard Kit: 200 Assorted Bands, 25 Sizes

Master Kit: 400 Assorted Bands, 25 Sizes

Reference Number

KIT15-050-21

KIT15-100-21

KIT15-200-21

KIT15-400-21



Lower Second Bicuspid, Plain

Starter Kit: 50 Assorted Bands, 27 Sizes

Intro Kit: 100 Assorted Bands, 27 Sizes

Standard Kit: 200 Assorted Bands, 27 Sizes

Master Kit: 400 Assorted Bands, 27 Sizes

Reference Number

KIT15-050-22

KIT15-100-22

KIT15-200-22

KIT15-400-22



Bicuspid Bands with Attachments

All GAC brackets and attachments are available pre-welded to our bicuspid bands.

Bands

Snap-Fit Cuspid Bands

Distinctive anatomical shaping characterizes GAC's line of cuspid bands in a total of 68 sizes: 19 upper, 15 lower, right and left. All Snap-Fit Bands are marked using our patented Electro Chemical Etching process for easy-to-read permanent identification, and the interiors are acid-etched for roughness to increase bond strength.



Snap-Fit Cuspid Band Kit Distribution - (UR, UL)

Band Size	000	OO	0	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16
Standard	4	4	4	4	4	4	6	8	8	8	8	8	6	4	4	4	4	4	4

Snap-Fit Cuspid Band Kit Distribution - (LR, LL)

Band Size	000	OO	0	1	2	3	4	5	6	7	8	9	10	11	12
Standard	4	4	6	6	6	8	8	8	8	8	8	8	6	6	6



Cuspid, plain

Standard Kit (200 assorted bands)

Reference Number

KIT14-200-30

Cuspid with Attachments

Any kit may be ordered with pre-welded attachments found in this catalog.

Snap-Fit Anterior Bands

GAC's anterior bands are precision-made to anatomical contours in a selection of 70 upper and lower sizes for easy placement and exact fit. In the maxillary there are 12 lateral and 16 central sizes, right and left. There are 14 sizes for the universally designed mandibular anteriors. All are made from seamless, stainless steel and polished to a smooth finish. All Snap-Fit Bands are marked using our patented Electro Chemical Etching process for easy-to-read permanent identification, and the interiors are acid-etched for roughness to increase bond strength.



Anterior, plain

Standard Kit (200 assorted bands)

Lower Only Intro Kit (100 assorted bands)

Reference Number

KIT14-200-00

KIT14-100-21

Anterior with Attachments

Any kit may be ordered with pre-welded attachments found in this catalog.

Band Material & Instruments

Band Material


Our two meter spools of stainless steel can be cut into strips and then pinched and welded to size.

Size	Reference Number
.10 x 4.6	26-321-06
.13 x 4.6	26-321-07
.15 x 4.6	26-321-08



Autoclavable Plastic Band Instruments

Instrument	Reference Number
Band Contouring Instrument	17-777-55
Bite Sticks: Contoured Point	17-777-44
Up Point	17-777-66
Down Point	17-777-77
Round Point	17-777-33
High Heat Bite Stick (Dry heat compatible)	17-777-22



Contouring Instrument



Bite Sticks
Contoured Point Up Point Down Point Round Point



High Heat Bite Stick

