



Extra Oral

DENTSPLY GAC offers a full line of facebows and accessories to meet your specific treatment needs. With a variety of choices for this particular segment, you will find popular facebows are available in a wide range of sizes, with or without cuspid

hooks and with or without white plastic coating. These can be paired with a selection of Class III traction and breakaway devices designed to provide safety and results for the clinician, and comfort and aesthetics for the patient.

For better dentistry

DENTSPLY
GAC


800.645.5530 | www.gacintl.com

Extra Oral Section Content

Facebows	241
Facebow Accessories	244
Class III Masks	245
Breakaway Traction Devices	246

Facebows

GAC offers a full line of facebows and accessories to meet your specific treatment needs. Standard facebows are available in a wide range of sizes, with or without cuspid hooks and with or without white plastic coating. All have a spring tempered outer bow securely welded to an inner .045" or .051" bow, and the size is clearly marked. Also available are Hamill Safety Facebows that have rounded distal ends to protect against intra and extra oral injury, as well as the Asher Lasher Facebow that is recommended by Dr. Cetlin and Roth as the facebow of choice for intrusion and retraction. CAUTION: Facebows should be used only with safety head/neck gear.



Roth .051" Facebows

Size	Reference Number	
	Special Short	Short
Universal	08-510-00	08-520-00
1	08-510-10	08-520-10
2	08-510-20	08-520-20
3	08-510-30	08-520-30
4	08-510-40	08-520-40
5	08-510-50	08-520-50
6	08-510-60	08-520-60

Facebows



Hamill Safety

The sharp distal ends of conventional facebow inner arches have caused numerous injuries, some severe. The rounded distal ends of our Hamill Safety Facebows protect against both intra and extra oral injuries. In addition, the easy spring action of the Hamill reduces the possibility of loose bands and makes insertion simple. Available in six sizes of both short and special short. Available only from GAC.

Developed by Dr. Maurice Hamill, Jr.

Inner Bow Size	Reference Number	
	Special Short	Short
1	08-419-10	08-429-10
2	08-419-20	08-429-20
3	08-419-30	08-429-30
4	08-419-40	08-429-40
5	08-419-50	08-429-50
6	08-419-60	08-429-60



Kahn Spur Hook Facebows

A special short facebow with a short spur centered on the inner bow.

Size	Reference Number
	Short
Universal	08-427-00
1	08-427-10
2	08-427-20
3	08-427-30
4	08-427-40
5	08-427-50
6	08-427-60



Asher Lasher Facebow

Recommended by Dr. Cetlin and Roth, the Asher Lasher Facebow is the facebow of choice for intrusion and retraction. With a high-pull cap, it retracts and intrudes; with a cervical strap, it simply retracts. The bow is hooked to the archwire mesial to the cuspid against a keyhole loop or a hook attached to the wire.

Description	Reference Number	
	Maxillary	Mandibular
Asher Lasher Facebow	03-100-99	03-100-89

Non Coated Facebows

Size	Reference Number				
	Special Short	Short	Medium	Long	High-pull
Universal	08-410-00	08-420-00	08-430-00	08-440-00	08-450-00
Universal/hook	08-412-00	08-422-00		08-442-00	
1	08-410-10	08-420-10	08-430-10		08-450-10
1/hook		08-422-10		08-442-10	
1.5	08-410-15	08-420-15			
2	08-410-20	08-420-20	08-430-20	08-440-20	08-450-20
2/hook	08-412-20	08-422-20		08-442-20	
3	08-410-30	08-420-30	08-430-30	08-440-30	08-450-30
3/hook	08-412-30	08-422-30			
3.5	08-410-35	08-420-35		08-440-35	
3.5/hook		08-422-35			
4	08-410-40	08-420-40	08-430-40	08-440-40	08-450-40
4/hook	08-412-40	08-422-40		08-442-40	
4.5	08-410-45	08-420-45		08-440-45	
4.5/hook		08-422-45			
5	08-410-50	08-420-50	08-430-50	08-440-50	08-450-50
5/hook	08-412-50	08-422-50		08-442-50	
5.5	08-410-55	08-420-55			
5.5/hook		08-422-55			
6	08-410-60	08-420-60	08-430-60	08-440-60	08-450-60
6/hook		08-422-60		08-442-60	
7	08-410-70	08-420-70	08-430-70		
8	08-410-80	08-420-80		08-440-80	



Plastic Coated Facebows

GAC offers a wide selection of facebows with an attractive white plastic coating. Coated facebows have proven to aid both patient comfort and cooperation, and can reduce irritation at the labial contact.

Size	Reference Number				
	Special Short	Short	Medium	Long	High-pull
Universal	08-413-00	08-423-00	08-433-00	08-443-00	
1		08-423-10			
1/hook	08-411-10				
2	08-413-20	08-423-20			08-453-20
3	08-413-30	08-423-30	08-433-30	08-443-30	08-453-30
3/hook	08-411-30				
3.5		08-423-35		08-443-35	
4	08-413-40	08-423-40	08-433-40	08-443-40	08-453-40
4/hook	08-411-40				
4.5		08-423-45			
5	08-413-50	08-423-50		08-443-50	08-453-50
5/hook	08-411-50				
5.5		08-423-55			
6	08-413-60	08-423-60			
8			08-433-80		



Facebow Accessories

Expansion Facebow

Distalizes and uprights molars simultaneously.



Size	Reference Number
Short	08-100-02

Micro Progressive Universal Facebow



Size	Reference Number
Short	08-700-01

Headgear Pouch (packs of 5)



	Reference Number	
Yellow	Pink	Teal Green
23-000-82	23-000-83	23-000-84

J-Hooks



	Reference Number	
	Small	Large
Pair	08-200-02	08-200-01

Facebow Stops



	Reference Number	
	.045"	.051"
Pack of two	08-229-21	08-229-22

Class III Masks

Protraction Face Mask

Centrally located to facilitate sleeping, our Protraction Face Mask is designed for high adjustability and patient comfort. Replacement parts are available for heavy wear, and the soft tissue contacts breathe for better patient hygiene.



Description	Reference Number
Protraction Face Mask	17-100-20
Replacement Pad and Key Kit	17-100-30

Reverse-Pull Headgear

GAC offers a universal Reverse-Pull Headgear with an elastic traction band. The built-in headband supports the mask and is made from soft, wide, elastic material for extra comfort. The forehead pad and chin cup are durable, adjustable, and non-irritating. All parts are easily adjustable to ensure an ideal fit. Another GAC exclusive.

Description	Reference Number
Mask, complete	17-100-00
Chin Cup	17-200-00
Head Pad	17-300-00
Elastic Strap	17-100-01
Face Bar with Screws	17-100-02



Ideal® Face Mask

Our zygoma-anchored Reverse-Pull Headgear is the ideal design for patient cooperation. The built-in traction band supports the mask. Every juncture is adjustable, and our padding is hypo-allergenic. The frame is highly stable, and all pads are designed to follow the anatomy. No pressure is placed on the chin.

Description	Reference Number
Ideal Face Mask	17-100-99



Breakaway Traction Devices

GAC offers a variety of breakaway traction devices, all designed to provide safety and results for the clinician and comfort and aesthetics for the patient.

SOF-Gear

Description

Reference Number

110 High-pull Cap	22-110-00
110 High-pull Cap - short	22-110-99
111 Chin Cup (for use with 110 caps)	22-111-00

Cervical Pad- 5" long	22-050-00
Cervical Pad- 7.5" long	22-075-00

Breakaway Cervical Set	23-100-00
Breakaway Velcro Modules	23-200-00
Velcro Hooks	22-000-12

UltraSoft™

Description

Reference Number

High-pull, Denim	KIT23-999-03
High-pull, Denim, no Modules	23-050-05
High-pull, Sunburst w/Modules	KIT23-000-03
High-pull, Sunburst, no Modules	23-050-07
High-pull, Denim w/Metal Modules	KIT23-888-03
High-pull, Sunburst w/Metal Modules	KIT23-111-03
High-pull, Starburst, Breakaway, Blue Metal Modules	KIT23-888-06
High-pull, Denim, Breakaway, Blue Metal Modules	KIT23-888-09

Cervical Headgear, Blue/Yellow	KIT23-999-01
Cervical Headgear, Red/Blue	KIT23-999-02
Cervical Neck Pad, Black	23-050-02
Cervical, Yellow, w/Metal Modules	KIT23-888-01
Cervical, Red, w/Metal Modules	KIT23-888-02
Cervical Headgear, Black w/Modules	KIT23-888-04

Breakaway, Yellow, Cervical, Blue Metal Modules	KIT23-888-07
Breakaway, Red, Cervical, Blue Metal Modules	KIT23-888-08
Breakaway Safety Modules, Plastic, pr.	23-050-06
Breakaway Safety Modules, Metal, Blue, pr.	23-050-08
Breakaway Safety Modules, Metal, Black, pr.	23-050-09
UltraSoft Neckpad, Red/Blue	23-050-03
UltraSoft Neckpad, Yellow/Blue	23-050-04
UltraSoft Neckpad, Denim	23-050-12
Safety Modules, pack of 2, Metal, Blue	23-050-11



NewGear™ Headgear

The NewGear from GAC is not only functional with its single latch snap, comfortable gel cushion, and precision force module, its sleek styling also provides patients with an attractive alternative to conventional headgear. The more appealing the appliance, the greater the likelihood of patient compliance and achieving the results you and your patients want. Available in red and black, cervical and high-pull.

Style	Reference Number	
	Red	Black
High-pull	23-400-03	23-400-02
Cervical	23-300-03	23-300-02

Stargear Headgear

The Stargear High-pull Headgear removes obstacles to patient cooperation with a design that does not interfere with a patient's hair style. The ponytail fits easily through the opening to allow the headgear to provide the comfort desired by the patient and the effective retraction required by the clinician. The kit includes the black Stargear High-pull cap with three sizing loops on each side and two black breakaway safety modules.

Description	Reference Number
Stargear High-pull Kit, black	KIT23-777-09
Stargear High-pull Cap, black	23-050-10
Stargear Cervical Kit, print	KIT23-777-01
Stargear Neckpad, print	23-050-01

