

Supplies

There are many choices when it comes to purchasing your orthodontic supplies. DENTSPLY GAC makes that choice simple, with an extensive selection to meet every need your practice may have. You'll find accessories to manage and organize your

bracket inventory, a wide array of bonding accessories, burs and patient care and safety items. Our goal is to help your office stay efficient by providing a complete line of quality products and supplies all from one source.

For better dentistry

DENTSPLY
GAC

800.645.5530 | www.gacintl.com

Supplies Section Content

Bracket Accessories	321
Bonding & Adhesive Accessories	324
Bur Accessories	330
Brace Relief Aids	333
Oral Hygiene	337
Retainers & Appliance Cases	340
Film & X-Ray Accessories	343

Bracket Accessories

Bracket Management System Base & Tower

Designed to hold, organize and store Mystique® MB, In-Ovation® R and In-Ovation® C boxes. Ergonomically designed system greatly facilitates bracket organization, bracket loading and inventory management.

System Benefits:

- Assistant friendly to reduce bracket setup time
- Saves valuable inventory management time
- Easy to read labeling
- Clear Bracket Tray Lid provides effortless inventory checks
- Reduces cross contamination concerns often associated with other delivery systems
- High grade material will stand up to most commonly used disinfectants and cleansers
- Convenient funnel for precise filling
- Tabletop presentation or wall mountable



	Reference Number
Base & Tower	20-300-00
Tray	20-300-01
System Tray Funnel	20-300-02
Credit Card Presentation	20-010-05
Bracket Management System Kit (1 each Base, Tower, Tray & Funnel)	KIT20-300-00

Credit Card Presentation Cards

Custom prescription set-ups can be ordered in this convenient credit card style presentation.



	Reference Number
Credit Card Presentation	20-010-05

Bracket Placement Cards

- Individual bracket trays with adhesive backing for temporary storage and setup
- Deep wells prevent bracket spillage
- Space is provided for the patients name, for easy identification



	Reference Number
Bracket Placement Card, Blue	20-050-99
Bracket Placement Card, Cardboard	20-050-00

Bracket Accessories

Bracket Tray Holder

- Manufactured from strong acrylic
- Holds up to 10 individual or disposable bracket trays
- Keeps setups away from contaminants and light
- Excellent for carrying bracket cases to satellite offices
- Attractive enough for any office décor



	Reference Number
Tray Holder	20-010-16

Consultation Models

- A complete line of orthodontic models for case presentations
- Allow for realistic examples of orthodontic appliance treatments



	Reference Number
Allure MB Consultation Model	TYP01-777-77
Balance Consultation Model	TYP31-777-77
Cetlin Consultation Model	TYP72-777-99
Elation MB Consultation Model	TYP133-777-77
In-Ovation C Consultation Model	TYP100-777-77
In-Ovation C/In-Ovation R Consultation Model	TYP100-777-89
In-Ovation R Alpern Consultation Model	TYP89-777-14
In-Ovation R Consult Model	TYP89-777-77
In-Ovation C Maxillary Model with Stand	TYP100-777-00
In-Ovation L Consultation Model	TYP90-777-77
In-Ovation MTM Consultation Model	TYP90-777-90
Micro Roth Consultation Model	TYP72-777-77
Mini Ovation Consultation Model	TYP85-777-77
Mini Shoulder Consultation Model	TYP75-777-77
Mystique MB/In-Ovation R Consultation Model	TYP00-777-89
Mystique MB Maxillary Model with Stand	TYP00-777-00
Mystique MB Consultation Model	TYP00-777-77
OmniArch Consultation Model	TYP77-777-77
OmniArch VS Consultation Model	TYP77-777-88
Ovation Consultation Model	TYP82-777-77
In-Ovation R Maxillary Consultation Model with Stand	TYP89-777-89
Single Wing Consultation Model	TYP72-777-88

Etchants

- Thin light cure gel
- Set time-unlimited (40 second light cure)
- Approximate number of bonds per kit – 400



	Reference Number
Ideal Conditioner 5mm Etch Syringe	52-777-03
Ideal Etching tips (pack of 8)	52-777-04

Mono Etch Compules

- Eliminates any chance of patient cross contamination during the etching process
- Each dispenser contains enough etch for 28 brackets
- The unique dispenser allows accurate etch placement on each tooth, eliminating waste
- Easy to handle, and use
- Contains 36% phosphoric acid
- Pack of 10 – 1.5g dispensers



	Reference Number
Mono Etch Compules	97-350-85

Prophy Angles and Paste

- Disposable angles are ideal for removing stains from occlusal pits and fissures and for preparing tooth surfaces prior to sealant and orthodontic applications
- Individual prophy paste cups for a germ-free environment



	Reference Number
Prophy Angles, disposable, pack of 144	97-300-50
Prophy Angle Cup Refills for 97-300-50 (144 cups)	97-350-80
Ideal Prophy Paste Packs, box of 200	97-155-55

Bonding & Adhesive Accessories

Absorb-Its™

- More effective for clinicians maintaining a dry field and more comfortable for the patient
- Flat, so they fit more easily into tight places
- Can be doubled-up for greater absorbency
- Made from the same super absorbent material used to make baby diapers
- This material is softer for greater patient comfort, and will not stick to the lining of the mouth
- Available in boxes of 800



	Reference Number
Absorb-Its™	97-111-12



Mixing Tools

- Intended for adhesives and cements which require preparation
- Ease of use and convenient package quantities

	Reference Number
Mixing Pads - 3" x 3" (pack of 10)	52-650-00
Mixing Pads - 3" x 3" (1 per pack)	51-000-07
Mixing Brushes pack of 50 - to apply sealants, etchants or primers	52-640-00
Mini sponge applicator - Wide bristles, thick handle - pack of 200	52-666-66

Ideal® Mini Mixing Sticks

Small design for easy adhesive mixing control and application



	Reference Number
Mini Mixing Stick	51-000-05

Micro Brushes Etch/Seal

- To apply sealants, etchants or primers; “sponge-like tip”
- Colors allow designation of tasks for prevention of cross contamination
- Pack of 100



	Reference Number
Purple	17-000-23
Green	17-000-24
Blue	17-000-25
Peach	17-000-26

Sterilizable Spatula

- Excellent for all Glass Ionomer cements
- Flexible for excellent control when mixing
- Rigid, easy handle grip prevents slippage
- Cleans easily, and will withstand repeated autoclaving and dry heat sterilization up to 370°F



	Reference Number
Spatula Sterilizable	17-555-55
Stainless Steel Cement Spatula No.24 Sm Blade	ORTHOPLI0241

Bite-ramps

- Works in open bite cases by allowing the anteriors to ride up and out of occlusion
- Swept design mimics the natural shape of the tooth - aids in the gentle redirection of anteriors
- High strength resin construction reduces trauma on enamel and increases patients comfort
- Pack of 10



	Reference Number
Bite-ramps	97-980-00

Bonding & Adhesive Accessories

Retract EEZ

The Retract EEZ is a clear, polycarbonate retractor that can be cold or heat sterilized. The angled end helps retract lips and is comfortable for the patient. The patient retracts lips with one hand and holds the mirror in the other, allowing the assistant to have both hands free to properly take photographs. One size fits all.



	Reference Number
Retract EEZ, each	97-999-10

Disposable Cheek Expanders

- Allows unrestricted access to the mouth by holding the mouth open and depressing the tongue to prevent injury
- Perfect for bonding, banding, and whitening procedures
- Manufactured from a flexible, easy to use plastic
- Available in three unique sizes to fit most patients
- Retractors are color coded by size for easy identification
- Available in packs of 10



	Reference Number
Small - Green	54-444-01
Medium - Blue	54-444-02
Large - Red	54-444-03

Open Wider Cheek Retractors

Autoclavable plastic Wider Cheek Retractors are available in regular size hard plastic and in regular and mini soft plastic.



	Reference Number
Regular, hard - clear	47-704-04
Regular, soft - light blue	47-704-05
Mini, soft - light blue	47-704-06

Barker Double Ended Retractor

- New large size available to fit most patients
- Each end of the retractor has a unique size
- The ends are elongated, resulting in increased access to the posterior region
- Pictures can be taken without the aid of an assistant, and the patient does not suffer any discomfort
- Available in cold sterilizable and autoclavable variations; 1 pair per pack



	Reference Number	Reference Number
Size	Cold Sterilizable	Autoclavable
Pedo Retractor	54-555-55	54-555-56
Large Retractor	54-222-02	54-222-01

BARKER is not a trademark of DENTSPLY GAC

Wyred™ Cheek Expander

- Spring wire support for lip and cheek retraction that provides consistent and adjustable force for improved access to the posterior teeth
- The bow is made from stainless steel, for repeated use
- Components can be ordered separately



	Reference Number
Retractors only - Heat Sterilizable 390° - Small	54-333-04
Retractors only - Heat Sterilizable 390° - Large	54-333-05

Disposable Tongue Guards

- Makes bonding easier for patients by relaxing the jaw during bonding procedures
- Keeps the tongue away from the teeth during bonding and banding procedures
- Soft plastic construction allows for easy trimming with scissors if needed
- Holes are provided to hold saliva ejectors for maintaining a dry field








	Reference Number
Red Cherry Scented - Small (pack of 50)	17-888-55
Red Cherry Scented - Large (pack of 50)	17-888-88

Bonding & Adhesive Accessories

Autoclavable Short Tipped Bitesticks

- Will withstand temperatures up to 260 degrees
- Down point stainless steel tip only
- Available in four bright colors, for work station organization



	Reference Number
 Autoclavable Bitestick - Teal	17-777-11
 Autoclavable Bitestick - Cobalt Blue	17-777-12
 Autoclavable Bitestick - Light Green	17-777-13
 Autoclavable Bitestick - Orange	17-777-14

Standard Tipped Heat Sterilizable Bitestick

- Will withstand temperatures up to 370 degrees



	Reference Number
Contoured Point	17-777-44
Up Point	17-777-66
Down Point	17-777-77
Round Point	17-777-33

Bracket Height Gauges

The GAC height gauge has permanent markings for precise placement of brackets from 3.5mm to 5mm from the incisal edge. Gauges are available in either .018 or .022 and they even hold brackets firmly for loading of adhesive. Heat sterilizable plastic or cold sterilizable aluminum.



	Reference Number	
	.018	.022
Aluminum	97-999-11	97-999-21
Plastic	97-999-12	97-999-22
Indirect bonding (aluminum) (universal)	97-999-15	

Disposable Mouth Mirrors

Approximately six inches in length and with a one inch diameter mirror, these disposable, white plastic mouth mirrors come with a peel-off, protective covering over the mirror. Lightweight, disposable, and economical.



	Reference Number
Disposable Mouth Mirrors, pack of 24	17-333-33

Photographic Mirrors

Photo mirrors are made from chromium - a hard, durable metallic element that can be autoclaved. These are available in adult size occlusal or in a kit of four mirrors: lingual, buccal, child occlusal, and adult occlusal.



	Reference Number
Photo Mirrors, set of 4	97-300-59
Photo Mirrors - Adult Occlusal	97-300-51

Carbide Fissure

- Cross-cut carbide fissure bur, safe- or standard-tipped available
- May be autoclaved, chemiclaved, or dry-heat sterilized - do not immerse in cold-sterilization solution
- Friction grip; Sold individually



Tip Type	Head Size	Head Length	Reference Number
Safe	0.9mm	5.2mm	25-699-55
Standard	0.9mm	5.2mm	25-699-99

Hilux 250 Curing Light

- This high quality light has all the features you are accustomed to paying more for: one touch trigger handpiece with built in timer, whisper quiet fan, versatile mounting ability, built in radiometer, and more
- Our exclusive Turbo Tip boosts output to up to 1200mw/cm², so curing is complete in just seconds
- Use with our Ideal® Light Cure CompuLe Kit for a complete bonding solution



	Reference Number
Hilux 250 Curing Light, 110 volt	97-500-10
Hilux Replacement Bulb, 12v/18w	97-500-01
Turbo Light Cure Tip	97-500-02
Hilux 250 Curing Light, 110 volt, w/Turbo Tip	KIT97-500-10

Hilux is not a trademark of DENTSPLY GAC

Bur Accessories

Debonding & Polishing Burs



Spiral-Fluted Debonding	25-140-02	25-140-03	25-140-01
Tip Type	Parallel Chamfer	Parallel Chamfer	Tapered Chamfer
Blade/Grit	12/Standard	12/Standard	12/Standard
Grip	Friction	Friction	Friction
Head Size	1.2mm	1.4mm	1.6mm
Head Length	8mm	8mm	8mm
Max Speed	<300,000rpm	<300,000rpm	<300,000rpm



Fluted Carbide 1.2mm Bur	25-140-20
Blade/Grit	20/Fine
Grip	Friction
Head Size	1.2mm
Head Length	8mm
Max Speed	<300,000rpm



Fluted Carbide Diamond Bur	25-780-29
Blade/Grit	12/Standard
Grip	Friction
Head Size	1.8mm
Max Speed	<300,000rpm



Universal Polisher

- Designed for all resin-based (composite) materials

	25-140-17
Angle	Right
Head Size	6mm
Head Length	10mm
Optimum Speed	5,000rpm

Diamond Ceramic Polisher

- Designed for all resin-based (composite) materials

	25-140-13	25-140-15
Description	Enamel Polish Disk	Knife edge wheel polisher
Angle	Right	Right
Gloss/Grit	High/Fine	High/Fine
Head Size	3.5mm	10mm
Head Length	10mm	2.5mm
Optimum Speed	5,000rpm-10,000rpm	5,000rpm-10,000rpm



Debonding/Enamel Polishing Kit

Includes:

- (3) Spiral-fluted debonding burs (Item # 25-140-01; 25-140-02 & 25-140-03)
- (1) Enamel Polishing Disc (Item # 25-140-11)
- (2) Diamond Impregnated Ceramic Polishers (Item # 25-140-13 & 25-140-15)
- (1) Universal Polisher (Item # 25-140-17)
- (1) Bur Block (Item # 25-141-90)



Description

Description	Reference Number
Debonding/Enamel Polishing Kit	25-140-00
Enamel Polishing Disc, 11mm	25-140-11

Bur Accessories

Magnetic Bur Block Autoclave

- Autoclavable Bur Block is heat sterilizable and holds 14 burs of varying sizes
- The block's magnetic property not only holds the burs securely within the block, but also holds the block in position on any metallic surface



	Reference Number
Magnetic Bur Block Autoclave	17-456-99

Cordless Dremel

- Cordless, straight dental lab hand dremel comes complete with charger and battery
- Special chuck included needed for dental /orthodontic burs and other attachments
- Two speeds - 5,000 RPM/10,000 RPM
- 90 day warrant against defective material or workmanship



	Reference Number
Cordless Dremel	25-990-01
Cordless Dremel Battery	25-990-02

Brace Relief Aids

Patient Relief Wax In Scented And Unscented


- Medical grade patient wax
- Ideal for temporary relief from appliance irritation
- The box design is small with rounded edges for patient convenience
- Available in scented and unscented color variations; Pack of 50 boxes



Description

Reference Number

Tropical Sparkle Blue	17-222-33
Tropical Sparkle Green	17-222-31
Tropical Sparkle Orange	17-222-34
Tropical Sparkle Pink	17-222-32
White/Unscented	17-222-22
Red Cinnamon/Scented	17-222-23



Silicone Patient Boxes

- Medical grade silicone used for temporary relief from appliance irritation
- Longer lasting and easier to use than patient wax
- Transparent and waterproof
- Resistant to hot liquid, and stays on the bracket longer for greater patient relief
- Pack of 50 boxes



Description

Reference Number

1 pack	17-222-00
Bulk pack of 75 silicone strips	17-222-05

Exerbite Therapy Wafer

- Reduces pain and discomfort
- Chewing on the exerbite wafer helps increase blood circulation around recently moved teeth reduce pain and discomfort
- Stimulates, massages and toughens the gum tissue
- Available in clear cherry; Pack of 25

Exerbite is not a trademark of DENTSPLY GAC



Description

Reference Number

Pack of 25	17-111-11
------------	-----------

Brace Relief Aids



Ortho-Lipguard™

- Provides temporary relief of direct pressure to the lips caused by Orthodontic appliances
- Applications include playing instruments and non contact sports
- Consists of a strip of medical grade plastic material that can be heated and formed to the contours of the patients mouth, creating a custom guard for each application
- The strip can be reheated and formed throughout orthodontic treatment, to compensate for moving teeth
- 5 strips per bag

Patent #6,080,923

Ortho-Lipguard is not a trademark of DENTSPLY GAC

Description

Ortho-Lipguard

Reference Number

17-222-24



Hook Cushions

- A great way to protect patients from unnecessary discomfort
- The cushion is placed over the bracket hooks, providing a buffer between the cheek and the bracket hook
- Tight fit eliminates slipping
- Soft, medical grade plastic does not swell or absorb water, making it safe and comfortable for the patient
- 14 cushions per stick, 10 sticks per pack - total of 140 cushions per package

Description

Hook Cushions

Reference Number

59-111-11



Breathable Brace Guards

- Innovative design provides upper and lower arch protection from braces while playing non contact sports
- Three breathing vents located at the front allow optimum airflow while playing sports
- Increased clearance for the tongue, which fits comfortably between the bite pads
- Available with or without a helmet strap
- Pack of 10

Description

Clear

Black

Green

Silver

Yellow

Multi-Color Pack

Reference Number

Without Strap

17-999-20

17-999-21

17-999-22

17-999-23

17-999-24

17-999-25

With Strap

17-999-30


17-999-31

17-999-32

17-999-33

17-999-34

17-999-35



Ideal® Brace Guards

- Dual Arch Mouthguard designed specifically for upper and lower arch protection
- Manufactured from a medical grade resin that can be heated and customized, making it soft enough to protect the gingival, but tough enough to stand up to hard use
- Available in vivid colors, allowing patients to match team colors
- Small and regular sizes available, with and without a strap; Packs of 10



Description

Regular

Description	Reference Number	
	With Strap	Without Strap
Yellow	17-999-00	17-999-03
Aqua	17-999-01	17-999-04
Silver	17-999-02	17-999-05
Clear	17-999-06	17-999-07

Small

Yellow	17-999-10	17-999-13
Aqua	17-999-11	17-999-14
Silver	17-999-12	17-999-15
Clear	17-999-16	17-999-17


Boil And Bite Mouthguards

- Allows for good protection and excellent air intake
- Includes a loop strap design that is easy to attach to helmet




Description

Description	Reference Number	
	Adult	Youth
Clear	25-198-96	25-198-86
Black	25-198-91	25-198-81
Blue	25-198-93	25-198-83
Red	25-198-94	25-198-84



Resin Mouthguard Sheets

- Resin Mouthguard Materials are used for children and young adults participating in low contact sports
- Variety of colors - coordinate with team colors!
- Box of 25 sheets



120"/3MM

Description	Reference Number	
	5" Square	125mm Circle
Clear	25-200-29	25-190-23

160"/4MM

Clear	25-200-30	25-190-24
Color	25-200-(Color #)	




Brace Relief Aids



Laminated Mouthguard Sheets

- Two layers laminated together equal higher tensile strength and uniform density
- 1/8" plastic insert protects incisal edges of teeth
- Latex Free
- 4mm Thickness
- Box of 12 sheets and 12 straps



Description

Reference Number

Black	25-200-04
Blue	25-200-03
Green	25-200-07
Orange	25-200-02
Purple	25-200-09
Red	25-200-05
White	25-200-06
Clear	25-200-01
Black/Orange	25-200-14
Black/Red	25-200-12
Blue/White	25-200-11
Blue/Yellow	25-200-10
Red/White	25-200-13
Red/White/Blue	25-200-16



Mouthguard Straps (Clear)

- Safety release allows the mouthguard to breakaway under 12 pounds of pressure
- Easily attaches to mouthguards using a Blazer Torch
- Pack of 25

Description

Reference Number

Mouthguard Straps	25-200-40
-------------------	-----------

Oral Hygiene

Retainer Retrievers

- Alleviates problems associated with removing clear appliances
- The patient can remove appliances from the buccal side cleanly and efficiently
- Small enough to store in an appliance case or on a key chain with the hole provided
- Packaged in individual tamper proof patient bags; 10 per pack

Patent #7,011,517




Description

Tropical Blue
Tropical Pink

Reference Number

97-000-00
97-000-01



Chewies™ Aligner Tray Series

- Simple solution for properly seating aligners with an imperfect fit
- Cleaner, more professional alternative to cotton rolls
- Provides uniform force in a specific area of the mouth where needed
- Available in packs of 10, tamper proof individual patient packs
- 2 chewies per patient pack




Description

Green/Mint
Pink/Bubble Gum
Purple/Grape
Yellow/Pineapple
White/Unscented

Reference Number

17-111-15
17-111-18
17-111-17
17-111-19
17-111-16



Aligner Patient Kit

The Aligner individual patient kit provides four key aligner products used during treatment, all in one handy kit. The products provided include the following:

- A package of two Chewies™ Aligner Tray Seaters – An essential product used to properly seat the tray after each appointment
- Two Retainer Retrievers – A great patient compliance aid that helps remove the tray from the arch quickly and efficiently
- One Econo Retainer Case – A durable retainer case that will fit almost any aligner appliance
- One Mini Retainer Case – A small case that's perfect for storing in a pocket or a purse
- All four products are essential for most aligner treatment
- Two color variations are available to cater to different patient requests



Description

Aligner Patient Kit Boys (Sparkle Cases: 1 Black Econo and 1 Green Mini, 2 White Retainer Retrievers and 2 Green Mint Aligner Tray Seaters)
Aligner Patient Kit Girls (Tropical Cases: 1 Pink Econo and 1 Blue Mini, 2 Green Retainer Retrievers and 2 Green Mint Aligner Tray Seaters)

Reference Number

97-950-56
97-950-55

Oral Hygiene



Retainer Brite® Cleaning Tablets

Save chair time when patients use Retainer Brite® Cleaning Tablets as they won't take up your time with appliance hygiene problems. Keep tablets on hand for in-office use! It's easy to supply your patients with Retainer Brite® Cleaning Tablets:

- Stock boxes of 1 or 3 month supplies in your office to sell or give to patients.
- Purchase Patient Paks and distribute to each patient when you deliver a new appliance. An easy order form is inside, so patients can purchase more Retainer Brite® Cleaning Tablets directly from DENTSPLY Raintree Essix.

Description

Reference Number

Retainer Brite® Cleaning Tablets 1 Month 36 Tablet Boxes	97-100-00
Retainer Brite® Cleaning Tablets Patient Pak (includes 2 tablets)	97-100-01
Retainer Brite® Cleaning Tablets 1 Month 36 Tablet Boxes (24 Boxes)	97-100-02
Retainer Brite® Cleaning Tablets Pack of 100	97-100-03



Interproximal Brushes

- Ideal for orthodontic patients
- Easy to use - small bristles easily clear around ortho appliances
- Individually packed for patient give-away
- Box of 100

Description

Reference Number

Interproximal Brushes	52-610-10
-----------------------	-----------



Proxi-Floss

Proxi-Floss is an interdental floss that will not snag, shred, or fray, and it is safe and easy for orthodontic patients to use. Available in bulk packs of 250 or a box of ten, 30 piece packets.

Description

Reference Number

Proxi-Floss, bulk pack of 250	25-793-00
Box of 10 packs, 30 pieces per pack	25-793-30

Ortho ER™

The Ortho ER emergency repair kit is designed to assist orthodontic patients in making repairs to their appliances, therefore avoiding unnecessary and disruptive office visits and after hours calls. The zippered ER case contains mouth and easel mirrors, travel toothbrush, wax, Soft Stick Wire Mender, and color instructions. Ortho ER is compact, easy to store, and inexpensive; and it is available in two, bright, translucent colors.



Description	Reference Number
Blue	97-950-61
Green	97-950-62

Description	Reference Number
Blue	97-950-61
Green	97-950-62

Remin+™ Gel


Remin+ is a multiple ion formulation that contains all the essential elements to harden teeth. Ideal for post tooth whitening, post interproximal stripping, and post occlusal equalibration. This mint gel is a patient take-home product that has the viscosity of a light toothpaste. Available in a 4oz squeeze bottle.



Description	Reference Number
Remin+ Remineralizing Solution, 4 oz bottle	25-745-00

Description	Reference Number
Remin+ Remineralizing Solution, 4 oz bottle	25-745-00

Retainers & Appliance Cases



Round Retainer Cases

- Unique round design is big enough to hold most retainers, but small enough to fit into a pocket or purse
- Superior locking mechanism that will stand up to heavy use

Description	Reference Number	
	1" - 10 Cases	1" - 300 Cases
Light Blue	12-000-01	KIT12-000-01
Navy Blue	12-000-02	KIT12-000-02
Green	12-000-03	KIT12-000-03
Hot Pink	12-000-04	KIT12-000-04
Red	12-000-05	KIT12-000-05
Black	12-000-06	KIT12-000-06



Sparkle Retainer Cases

- Sleek 1.2" x 3.4" modern design is big enough to hold most retainers, but small enough to fit into a pocket or purse
- Superior locking mechanism that will stand up to heavy use
- Available in packs of 10 cases

Description	Reference Number	
	10 Cases	330 Cases
Sparkle Blue	12-100-02	12-330-02
Sparkle Green	12-100-03	12-330-03
Sparkle Red	12-100-05	12-330-05
Sparkle Black	12-100-06	12-330-06
Sparkle Purple	12-100-07	12-330-07
Sparkle Gold	12-100-08	12-330-08
Sparkle Glow	12-100-09	12-330-09
Sparkle Plum	12-100-10	12-330-10
Sparkle Cool Mint/Aqua	12-100-11	12-330-11
Tropical Sparkle Blue	12-100-12	
Tropical Sparkle Green	12-100-13	
Tropical Sparkle Orange	12-100-14	
Tropical Sparkle Pink	12-100-15	
Glow Orange	12-100-21	
Glow Pink	12-100-22	
Glow Green	12-100-23	KIT12-100-23

Bionator Cases

- Extra deep - 1½" deep x 3" wide by 2½" long
- Strong hinge and secure lock protects appliances sufficiently
- Available in packs of 10 cases

Bionator is not a trademark of DENTSPLY GAC



Description

Sparkle Black Case

Reference Number

12-200-04

Lab Appliance Containers

- 2" deep sizes
- The top of each case is vented, reducing unpleasant odors
- Available in packs of 6 cases

Description

Navy Blue

Sparkle Green

Sparkle Orange

Sparkle Pink

Reference Number

12-200-02

12-200-05

12-200-06

12-200-07



GermGuard Retainer Cases

- Designed to kill bacteria that can spread retainer cases between regular cleaning
- Unique updated design

Description

Sparkle Black

Sparkle Blue

Reference Number








12-300-01

12-300-02



Retainers & Appliance Cases



-  Sparkle Blue
-  Sparkle Green
-  Sparkle Black
-  Tropical Sparkle Blue
-  Tropical Sparkle Green
-  Tropical Sparkle Orange
-  Tropical Sparkle Pink

Mini Retainer Cases

- Designed with sleek 1" high, 2" wide, and 3.375" long rounded edges for ultimate carrying comfort
- Great for aligners and thermoformed retainers

Description	Reference Number	
	25 Cases	500 Cases
Sparkle Blue	12-400-02	12-500-02
Sparkle Green	12-400-03	12-500-03
Sparkle Black	12-400-06	12-500-06
Tropical Sparkle Blue	12-400-12	12-500-12
Tropical Sparkle Green	12-400-13	12-500-13
Tropical Sparkle Orange	12-400-14	12-500-14
Tropical Sparkle Pink	12-400-15	12-500-15

Economy Retainer Case

- 1" pack of 300

Description	Reference Number
Blue	12-555-55
Retainer Case Labels 100	12-000-00



Film & X-Ray Accessories

Orthomat™ X-Ray Film

You have to have x-ray film, but you don't have to pay a lot for it. GAC's Orthomat x-ray film is high quality, high contrast, blue or green sensitive film that provides consistent results and sharp detail. Our x-ray film is packaged in 100 sheet packs and is available an affordable price.

Size	Reference Number		
	High Speed, Blue	Medium Speed, Green	
8" x 10"	17-545-00	17-545-99	
5" x 12"	17-380-00	17-380-99	
15cm x 30cm	17-390-00	17-390-99	
18cm x 24cm	17-350-00	17-350-99	
	800 Speed	400 Speed	100 Speed
Cassettes, 1 pr, 8" x 10"	17-391-00	17-391-99	17-391-55



Agfa Curix® Duplicating Film

Agfa Curix Duplicating Film provides especially accurate reproduction of all diagnostic images. It is available in both blue tinted and clear polyester bases, which provide excellent handling characteristics and exceptionally long shelf life. Available in boxes of 100 sheets, Agfa Curix Duplicating Film offers great, long lasting reproductions.

Description	Reference Number	
	Blue	Clear
Duplicating Film: 8" x 10"	17-370-00	17-370-05
5" x 12"	17-370-50	



Cephalometric Viewer Light Box

Our all metal Cephalometric Light Box offers a full 8.5" x 11" viewing surface with a standard daylight fluorescent lamp. It is small enough to be used on a desk or countertop or conveniently carried in a briefcase. Its leg bar provides easy ceph viewing, yet it swings away for low profile storage.

Description	Reference Number
Cephalometric Light Box	97-600-99



Film & X-Ray Accessories



Cephalometric Tracing Paper

This 100-page tablet contains 8" x 10" acetate sheets that resist tearing and curling. Each sheet is .003" thick and is protected by tissue paper inserts. Also available in kit form with four colored pencils and a sharpener.

Description

Reference Number

Cephalometric Tracing Paper Tablet

17-222-11

Kit: Includes Tracing Paper Tablet, 4 Pencils and Sharpener

17-222-10

Orthodontic Rulers

Description

Reference Number

Overbite Ruler

07-450-00

Orthodontic Ruler and Ceph Tracing

97-600-50

Archwire Computer

49-000-00

