

In-Ovation® R Offset Bicuspid Brackets

easy, fast, and convenient



Kate and Michael: actual In-Ovation patients.


Reduce office visits and speed up treatment time...

- Bonds easily to teeth
- Bracket only partially erupted teeth
- Improved comfort
- Easy to maintain oral hygiene

New In-Ovation R gingivally offset bicuspid brackets for correcting partially erupted bicuspids

Introducing In-Ovation's new offset lower bicuspid brackets. The pad's contour ideally fits the tooth surface and is positioned gingivally which reduces the likelihood of upper cusp interference. This means easy bonding and fewer office visits to repair broken brackets.

The offset design will also facilitate earlier ideal bracket placement on erupting bicuspids which may shorten overall treatment time. Easy for patients to keep clean in and around the gum line, oral hygiene is better maintained as compared to bands.



Easy to use

Precise anatomical fit for easy bonding as well as an easy to understand distogingival color code ID system.

Greater efficiency

Shorter treatment time as clinicians do not have to wait for slow erupting bicuspids to emerge. Fewer repair visits because of the bracket's offset design mean less chairtime and a more optimized practice.

Increased patient comfort

Less gingival irritation and easier to maintain good oral hygiene as compared with bands.

Product class: 89

Packaging: Packs of 10 (or included in patient KIT setup)

Item numbers:

89-34X-08 / Inov R ROT Lower Right 4 -17T-1A 4D OB

89-44X-08 / Inov R ROT Lower Left 4 -17T-1A 4D OB

89-35X-08 / Inov R ROT Lower Right 5 -22T-1A 4D OB

89-45X-08 / Inov R ROT Lower Left 5 -22T-1A 4D OB

For more information, call 800.645.5530.



The Science of Innovation. The Art of Patient Care.

www.gacintl.com