

Your Total Orthodontic Solutions Partner...

As an industry leader, GAC is your total orthodontic supplier offering a complete range of solutions that can easily and cost-effectively be incorporated into your existing practice and treatment plans. At GAC, we believe in leading innovation. That's why we provide state-of-the-art products from brackets, force-intuitive wires and highly efficient buccal tubes to bonding, auxiliary appliances and innovative software programs. Our mission is to help you achieve rapid, optimal treatment results while streamlining your orthodontic practice. Please ask your GAC representative for more information, call 800.645.5530, or visit www.gacintl.com.

Stylish twin systems of performance and strength.



OmniArch®



MicroArch®





© 2008 GAC International LLC. All rights reserved. 5.2008-300 - U.S. Patents: 4,415,530 and 4,889,485. Other U.S. and foreign patents pending and/or approved. 120-777-01

OmniArch®

The Most Advanced Versatile Twin Design

OmniArch is the first mid-sized bracket system that combines all the outstanding elements of the proven twin system with the most advanced bracket design available. Clinicians gain superior control, wedged torque in base, easy bonding; patients benefit from a bracket that is comfortable, aesthetically pleasing, and from outstanding results. Conveniently, OmniArch is available in a wide range of prescriptions and options.



System Benefits

- Aligned angulation for unequalled mesiodistal/incisogingival orientation
- Wedged torque for aligned slots and placement advantages
- Low profile provides greater comfort
- Mid-sized and rhomboid shaped for better performance and aesthetics
- GAC's SuperMesh® base for consistent bonding and fast clean-up


Prescription and Option Choices

- Roth Rx
- Roncone Rx
- Euro Rx
- Offset Bicuspids
- Andrews (7/3) Rx
- Golden OmniArch in a titanium nitride finish for aesthetics and biocompatibility
- MicroProgressive with high torque anteriors
- V-Slot in Cetlin® or BiDimensional Rx

MicroArch®

The Advanced Low-profile Twin

With many of the same features as OmniArch, MicroArch has proven to be one of the most popular mini brackets in the world for both clinicians and patients. A leading solution in performance, comfort and strength, each bracket is milled from aerospace-strength steel to avoid variations associated with casting and injection molding. To prevent electrolytic reactions, the base and bracket are brazed with a special gold alloy. To enhance patient comfort, the MicroArch has an exceptionally low profile and chamfered tie wings; while its SuperMesh® base means reliable bonding and debonding. The result is an aesthetic bracket that excels in both performance and reliability.



System Benefits

- Small size and rhomboid shape for greater comfort and aesthetics
- Easy archwire insertion
- Fast and reliable bonding with SuperMesh® base
- Smooth swept tie wings for greater patient comfort
- In/Outs built-in, not cut, for greater precision

Prescription and Option Choices

- Roth Rx
- Andrews (7/3) Rx
- MicroArch SB™
- Single Wing MicroArch
- Offset Base Bicuspids
- Standard Edgewise MicroArch
- Microprogressive Appliances
- V-Slot in Cetlin® or BiDimensional Rx

