

Tray Aligner Accessories

simple, convenient, easy to use



Solutions that support clinicians using tray aligners

- Ideal® Interproximal Reduction System
- Photo Mirrors
- Barker Photo Retractors
- Chewies™ Aligner Tray Seaters
- Retainer Retrievers
- Sparkle Retainer Cases
- Aligner Patient Kits

Tray Aligner Accessories: affordable solutions for clinicians and their patients

Ideal® Interproximal Reduction System

- Flexible strips
- Outstanding reliability
- Safe and precise



IDEAL is your one-stop interproximal stripping solution. With tray aligners, such as GAC's Cosmetic Correctors™, the judicious stripping of enamel is occasionally needed for effective treatment. IDEAL is designed to be gentle. A fine diamond coating on the first strip prevents catching or gouging during the first reduction phase, and facilitates backing out. Succeeding strips have varying grades of coating for all phases of interproximal stripping; exceptional flexibility contours along the tooth surface for the highest precision.

KIT04-100-00 Ideal® Interproximal Reduction System Kit

Photo Mirrors

Available in adult size occlusal or in a kit of four mirrors: lingual, buccal, child occlusal, and adult occlusal.

Item 97-300-51 Photo mirror, adult occlusal
Item 97-300-59 Set of four mirrors



Barker Photo Retractors

Comfortable, swept design allows the clinician or assistant to comfortably hold the patient's soft tissue away while using photography to document cases. Available in cold sterilizable or autoclavable. Two retractors per pack.



Item 54-555-55 Regular (Clear)
Item 54-555-56 Autoclavable (Gray)



Chewies™ Aligner Tray Seaters

- Simple solution for properly seating aligners with an imperfect fit
- Cleaner, more professional alternative to cotton rolls
- Provides uniform force in a specific area of the mouth
- Ten per pack in individual tamper proof patient packs
- 2 chewies per patient pack
- Available in five colors

Item 17-111-16 (White)




Retainer Retrievers

- Alleviates problems with removing clear appliances
- The patient can remove appliances from the buccal side cleanly and efficiently
- Small enough to store in an appliance case or on a key chain with the hole provided
- Ten per pack in individual tamper proof patient bags

Item 97-000-00 (Tropical Blue)


Item 97-000-01 (Tropical Pink)



Sparkle Retainer Cases

- Sleek 1.2" x 3.4" modern design is big enough to hold most retainers, but small enough to fit into a pocket or purse
- Superior locking mechanism that will stand up to heavy use
- Available in packs of 10 cases
- Available in 16 colors

Item 12-100-08 (Sparkle Gold)



Aligner Patient Kit

- Provides four key aligner products essential for use in most aligner treatment cases, all in one handy kit, including:
 - * A package of two Chewies™
 - * Two Retainer Retrievers
 - * One Econo Retainer Case
 - * One Mini Retainer Case

Item 97-950-55 (Pink/Blue/Green)

Item 97-950-56 (Black/Green/White)



The Science of Innovation. The Art of Patient Care.

www.gacintl.com