

Dr. \_\_\_\_\_  
 Address \_\_\_\_\_  
 City & State \_\_\_\_\_  
 Tel \_\_\_\_\_ Fax \_\_\_\_\_  
 Patient \_\_\_\_\_  
 Shipped Date \_\_\_\_\_ Placement Date \_\_\_\_\_  
 Date Received by Lab \_\_\_\_\_

PLEASE WRITE SPECIAL INSTRUCTIONS

\_\_\_\_\_

\_\_\_\_\_

\_\_\_\_\_

\_\_\_\_\_

\_\_\_\_\_

\_\_\_\_\_

\_\_\_\_\_

\_\_\_\_\_

PLEASE CALL ME ABOUT THIS ORDER

**GACOrthoLab**  
orthodontic laboratory services

**2525 3 MILE ROAD  
RACINE, WI 53404**

Toll-FREE: (866) 463-4300 Main Tel: (262) 752-4040  
 Fax: (262) 752-4060 For billing questions: (800) 645-5530

**Anterior Repositioning Splint**  
 Upper  Lower

Full coverage upper appliance with anterior ramp designed to "capture" the lower anterior teeth and hold the mandible in a forward/anterior position. Lower appliance consists of posterior coverage indexed with upper lingual cusps connected with an anterior lingual band of acrylic. Forward upper and lower models. Bite reflecting anterior position.

**Superior Repositioning Splint**  
 Upper (Roth Style)  Lower (Tanner splint)

Full coverage on upper arch with anterior ramp designed to provide protrusive and cuspid protection. Both upper and lower models with centric relation bite required. Bite should provide 1-2 mm minimum anterior opening for acrylic durability.

Mounted models for adjustable articulators are recommended.

**Flat Plane Occlusal Splint**  
 Upper  Lower

Full coverage with even contact of opposing arch. Counter model suggested with wax bite. Used as bruxism or Pre treatment diagnostic splint.

**MORA Splint**

Fabricated on lower arch with coverage of posterior teeth which includes upper posterior lingual cusp tip stops. Standard design includes heavy lingual connecting bar.

**Gelb Splint**

Identical to MORA splint, but includes anterior lingual band of acrylic.

**Materials**

Acrylic  Clear  Color \_\_\_\_\_

Pressure formed Biocryl base  Add Acrylic  No Acrylic

Soft Vinyl appliance - **Bruxism or Sports**

Dual Hardness - Soft inner surface with hard outer shell

**Splint Construction Bite**

Articulated Models Enclosed

- SAM
- Denar
- Panadent High Post
- Panadent Low Post
  - Quick Split Magnets  Upper only  Lower only
  - Other \_\_\_\_\_

Use Enclosed Wax Construction Bite

Articulate Models in Lab with Average Vertical Opening

**Clasps**

Mark location on diagram

Ball  Adams

Arrow  "C"

**PLEASE SHIP EXTRA:**

SHIPPING BOXES

PRE-PAID BAGS

PRESCRIPTION SHEETS

