


Extraoral Products

Facebows

GAC offers a full line of facebows and accessories to meet your specific treatment needs. Standard facebows are available in a wide range of sizes, with or without cuspid hooks and with or without white, plastic coating. All have a spring tempered outer bow securely welded to an inner .045" or .051" bow, and the size is clearly marked. Also available are Hamill Safety Facebows that have rounded distal ends to protect against intra and extra oral injury, as well as the Asher Lasher Facebow that is recommended by Dr. Cetlin and Roth as the facebow of choice for intrusion and retraction. **CAUTION: Facebows should be used only with safety head/neck gear.**



Roth .051" Facebows		
Size	Reference Number	
	Special Short	Short
Universal	08-510-00	08-520-00
1	08-510-10	08-520-10
2	08-510-20	08-520-20
3	08-510-30	08-520-30
4	08-510-40	08-520-40
5	08-510-50	08-520-50
6	08-510-60	08-520-60

Size	Reference Number							
	Special Short		Short		Medium	Long		High Pull
	Plain	Plain w/hk	Plain	Plain w/hk	Plain	Plain	Plain w/hk	Plain
1	08-410-10		08-420-10	08-422-10	08-430-10		08-442-10	08-450-10
1.5	08-410-15		08-420-15					
2	08-410-20	08-412-20	08-420-20	08-422-20	08-430-20	08-440-20	08-442-20	08-450-20
3	08-410-30	08-412-30	08-420-30	08-422-30	08-430-30	08-440-30		08-450-30
3.5	08-410-35		08-420-35	08-422-35		08-440-35		
4	08-410-40	08-412-40	08-420-40	08-422-40	08-430-40	08-440-40	08-442-40	08-450-40
4.5	08-410-45		08-420-45	08-422-45		08-440-45		
5	08-410-50	08-412-50	08-420-50	08-422-50	08-430-50	08-440-50	08-442-50	08-450-50
5.5	08-410-55		08-420-55	08-422-55				
6	08-410-60		08-420-60	08-422-60	08-430-60	08-440-60	08-442-60	08-450-60
7	08-410-70		08-420-70		08-430-70			
8	08-410-80		08-420-80			08-440-80		
UNI	08-410-00	08-412-00	08-420-00	08-422-00	08-430-00	08-440-00	08-442-00	08-450-00

Extraoral Products



Facebows

Hamill Safety

The sharp distal ends of conventional facebow inner arches have caused numerous injuries, some severe. The rounded distal ends of our Hamill Safety Facebows protect against both intra and extra oral injuries. In addition, the easy spring action of the Hamill reduces the possibility of loose bands and makes insertion simple. Available in six sizes of both short and special short. Available only from GAC.

Developed by Dr. Maurice Hamill, Jr.

Inner Bow Size	Reference Number	
	Special Short	Short
1	08-419-10	08-429-10
2	08-419-20	08-429-20
3	08-419-30	08-429-30
4	08-419-40	08-429-40
5	08-419-50	08-429-50
6	08-419-60	08-429-60



Kahn Spur Hook Facebows

A special short facebow with a short spur centered on the inner bow.

Size	Reference Number	
	Universal	Short
1	08-427-10	08-427-10
2	08-427-20	08-427-20
3	08-427-30	08-427-30
4	08-427-40	08-427-40
5	08-427-50	08-427-50
6	08-427-60	08-427-60



Asher Lasher Facebow

Recommended by Dr. Cetlin and Roth, the Asher Lasher Facebow is the facebow of choice for intrusion and retraction. With a high-pull cap, it retracts and intrudes; with a cervical strap, it simply retracts. The bow is hooked to the arch wire mesial to the cuspid against a keyhole loop or a hook attached to the wire.

Asher Lasher Facebow	Reference Number	
	Maxillary	Mandibular
	03-100-99	03-100-89



Facebow & Accessories

Plastic Coated Facebows

GAC offers a wide selection of facebows with an attractive white, plastic coating. Coated facebows have proven to aid both patient comfort and cooperation, and can reduce irritation at the labial contact.


Size	Reference Number					
	Special Short	Special Short w/cuspid hk	Short	Medium	Long	High Pull
Univ	08-413-00		08-423-00	08-433-00	08-443-00	
1		08-411-10	08-423-10			
2	08-413-20		08-423-20			08-453-20
3	08-413-30	08-411-30	08-423-30	08-433-30	08-443-30	08-453-30
3.5			08-423-35		08-443-35	
4	08-413-40	08-411-40	08-423-40	08-433-40	08-443-40	08-453-40
4.5			08-423-45			
5	08-413-50	08-411-50	08-423-50		08-443-50	08-453-50
5.5			08-423-55			
6	08-413-60		08-423-60			
7						
8				08-433-80		

Facebow Accessories

Expansion Facebow

Distalizes and uprights molars simultaneously.

Size	Reference Number
Short	08-100-02



Micro Progressive Universal Facebow

Size	Reference Number
Short	08-700-01



Headgear Pouch (packs of 5)

	Yellow	Pink	Teal Green
Reference Number	23-000-82	23-000-83	23-000-84

J-Hooks

	small	large
Reference Number	08-200-02	08-200-01



Facebow Stops

	.045"	.051"
Reference Number	08-229-21	08-229-22

Package of two



Extraoral Products

Class III Masks



Protraction Face Mask

Centrally located to facilitate sleeping, our Protraction Face Mask is designed for high adjustability and patient comfort. Replacement parts are available for heavy wear, and the soft tissue contacts breathe for better patient hygiene.

Description	Reference Number
Protraction Face Mask	17-100-20
Replacement Head Pad	17-100-21
Replacement Chin Cup	17-100-22



Reverse-Pull Headgear

GAC offers a universal Reverse-Pull Headgear with an elastic traction band. The built-in headband supports the mask and is made from soft, wide, elastic material for extra comfort. The forehead pad and chin cup are durable, adjustable, and non-irritating. All parts are easily adjustable to ensure an ideal fit. Another GAC exclusive.

Description	Reference Number
Mask, complete	17-100-00
Chin Cup	17-200-00
Head Pad	17-300-00
Elastic Strap	17-100-01
Face Bar with Screws	17-100-02



Ideal® Face Mask

Our zygoma-anchored reverse-pull headgear is the ideal design for patient cooperation. The built-in traction band supports the mask. Every juncture is adjustable, and our padding is hypo-allergenic. The frame is highly stable, and all pads are designed to follow the anatomy. No pressure is placed on the chin.

Description	Reference Number
Ideal Face Mask	17-100-99

Breakaway Traction Devices

Breakaway Traction Devices

GAC offers a variety of breakaway traction devices, all designed to provide safety and results for the clinician and comfort and aesthetics for the patient.

SOF-Gear

Description	Reference Number
Breakaway Cervical Set	23-100-00
Breakaway Velcro Modules	23-200-00
Velcro Hooks	22-000-12
110 High Pull Cap	22-110-00
110 High Pull Cap - short	22-110-99
111 Chin Cup (for use with 110 caps)	22-111-00
Cervical Pad- 5" long	22-050-00
Cervical Pad- 7.5" long	22-075-00



22-110-00 22-111-00

UltraSoft™

Description	Reference Number
Cervical Headgear, Blue/Yellow	KIT23-999-01
Cervical Headgear, Red/Blue	KIT23-999-02
Hi Pull - Headgear, Denim	KIT23-999-03
Hi Pull - Denim, no Modules	23-050-05
Hi Pull - Sunburst w/Modules	KIT23-000-03
Hi Pull - Sunburst, no Modules	23-050-07
Breakaway Safety Modules, Plastic, pr.	23-050-06
Breakaway Safety Modules, Metal, Blue, pr.	23-050-08
Breakaway Safety Modules, Metal, Black, pr.	23-050-09
Cervical Neck Pad, Black	23-050-02
Cervical, Yellow, w/Metal Modules	KIT23-888-01
Cervical, Red, w/Metal Modules	KIT23-888-02
High Pull Denim w/Metal Modules	KIT23-888-03
Cervical Headgear, Black w/Modules	KIT23-888-04
High Pull, Sunburst w/Metal Modules	KIT23-111-03
Breakaway, Yellow, Cervical, Blue Metal Modules	KIT23-888-07
Breakaway, Red, Cervical, Blue Metal Modules	KIT23-888-08
High Pull, Starburst, Breakaway, Blue Metal Modules	KIT23-888-06
High Pull, Denim, Breakaway, Blue Metal Modules	KIT23-888-09
Safety Modules, pk/2 Metal, Blue	23-050-11
UltraSoft Neckpad, Red/Blue	23-050-03
UltraSoft Neckpad, Yellow/Blue	23-050-04
UltraSoft Neckpad, Denim	23-050-12




23-100-00



KIT23-999-01



KIT23-999-03




23-050-12

23-050-11

23-000-84

Extraoral Products

Breakaway Traction Devices



NewGear™ Headgear

The NewGear from GAC is not only functional with its single latch snap, comfortable gel cushion, and precision force module, its sleek styling also provides patients with an attractive alternative to conventional headgear. The more appealing the appliance, the greater the likelihood of patient compliance and achieving the results you and your patients want. Available in red and black, cervical and high pull.

Style	Reference Number	
	Red	Black
Cervical	23-300-03	23-300-02
High Pull	23-400-03	23-400-02



Stargear Headgear

The Stargear Highpull Headgear removes obstacles to patient cooperation with a design that does not interfere with a patient's hair style. The ponytail fits easily through the opening to allow the headgear to provide the comfort desired by the patient and the effective retraction required by the clinician. The kit includes the black Stargear Highpull cap with three sizing loops on each side and two black breakaway safety modules.

Description	Reference Number
Stargear Highpull Kit, black	KIT23-777-09
Stargear Cervical Kit, print	KIT23-777-01
Stargear Neckpad, print	23-050-01
Stargear Highpull Cap, black	23-050-10